



# 2020 CATALOG

Dental and Veterinary Dental Products

---

IMPROVING ORAL HEALTH THROUGH  
ANTIMICROBIAL SELDOX® TECHNOLOGY





## DenteShield®

### PIT AND FISSURE SEALANT

DenteShield® Pit and Fissure Sealant works harder than traditional sealants by, actively killing bacteria on contact and preventing plaque and decay both on and beneath the sealant. DenteShield® Pit and Fissure Sealant PLUS adds fluoride, for plaque protection PLUS enamel repair.



DenteShield® Pit & Fissure Sealant PLUS - Clear (25 pack)

SKU: SBD-650-3



DenteShield® Pit and Fissure Sealant Clear Kit

SKU: SBD-600-1



DenteShield® Pit and Fissure Sealant Opaque Kit

SKU: SBD-600-5

## DenteShield®

### PRIMER

DenteShield® Primer is the first FDA-cleared primer to offer both a physical barrier to bacteria and acid as well as SELDOX® antimicrobial technology to significantly reduce the occurrence of white spot lesions. A bond-enhancing light cured resin, DenteShield® Primer is applied at the start of treatment directly after etching and works with your current adhesive.



DenteShield® Primer Light Cure - 3ml

SKU: SBD-200-5



DenteShield® Primer Light Cure - 7ml

SKU: SBD-200-1

## DenteShield®

### SURFACE SEALANT

DenteShield® Surface Sealant is a general use, light cured, unfilled resin containing SELDOX® antimicrobial technology to protect enamel or dentin non-occlusal surfaces against bacteria and plaque formation. And for orthodontists who prefer not to change bonding techniques. Surface Sealant may be applied to protect surface areas around brackets at the start of or mid-treatment.



DenteShield® Surface Sealant Light Cure - 3ml

SKU: SBD-100-5



DenteShield® Surface Sealant Light Cure - 7ml

SKU: SBD-100-1



GingiShield® is a clear, light-cured resin sealant applied by a veterinary professional at the end of a cleaning and can last a full year. The sealant, with antimicrobial SELDOX®, neutralizes the bacteria responsible for plaque formation and the development of periodontal disease.



## GingiShield® SEALANT STARTER KIT

The GingiShield® Antibacterial Enamel Sealant Starter Kit includes everything needed to begin protecting your patient's teeth with GingiShield®.

Designed for veterinary practices that do not currently have a curing light, this kit includes:

- High Speed Cordless LED Curing Light
- 2 x 1.5ml Sealant Syringes
- 2 x 1.2ml Etchant Syringes (37% phosphoric acid)
- 100ml Fluoride-free Prophy Paste
- Programmable Digital Timer for Etching
- Syringe Tips, Mixing Wells, Applicator Brushes
- Instructions For Use
- 30 Tri-Fold Client Brochures

GingiShield® Sealant Starter Kit  
SKU: SBD-950-2



## GingiShield® SEALANT KIT

The GingiShield® Antibacterial Enamel Sealant Kit is for veterinary professionals already using light-cured dental materials, or as a refill kit.

The GingiShield® Antibacterial Sealant Kit includes:

- 2 x 1.5ml Sealant Syringes
- 2 x 1.2ml Etchant Syringes (37% phosphoric acid)
- Syringe Tips, Mixing Wells, Applicator Brushes
- Instructions For Use

*Each kit will provide roughly 20 patient treatments.*

GingiShield® Sealant Kit  
SKU: SBD-950-6



## SelenBio Dental

SelenBio Dental, Inc. is a biotechnology company specializing in products that prevent the formation of dental biofilms.

### SELDOX® Antimicrobial Technology

SELDOX® is an antimicrobial technology able to inhibit bacterial biofilm formation on treated surfaces. SELDOX® is covalently linked, catalytic, non-toxic, and does not leach into the mouth or surrounding environment. These properties allow for effective and safe bacteria and plaque prevention, and improved oral health.

SelenBio Dental's sealant materials are shown to aid in the prevention of:

- Plaque
- Dental caries (cavities)
- Demineralization and white spot lesions
- Recurrent decay
- Restoration failure
- Gingival inflammation

*Products are FDA-cleared and have undergone successful university clinical trials.*



4425 Mopac Expwy South  
Building II, Suite 300A  
Austin, Texas 78735

512.887.2050

[info@selenbiodental.com](mailto:info@selenbiodental.com)